

OpenSees: Static, Cyclic, and Dynamic Analysis Examples

Frank McKenna

OpenSees Users Workshop

Berkeley, CA

August 26, 2005


Sponsored by the National Science Foundation
through the Pacific Earthquake Engineering Research Center
and through NEESit



Model (example.tcl):

```
model Basic -ndm 2 -ndf 3
node 1 0.0 0.0
node 2 504.0 0.0
node 3 0.0 432.0
node 4 504.0 432.0
fix 1 1 1 1
fix 2 1 1 1
mass 3 5.18 0.0 0.0
mass 4 5.18 0.0 0.0
geomTransf Linear 1
element elasticBeamColumn 1 1 3 3600.0 4227.0 1080000.0 1
element elasticBeamColumn 1 4 3 3600.0 4227.0 1080000.0 1
element elasticBeamColumn 3 3 4 5760.0 4227.0 4423680.0 1
pattern Plain 1 Linear {
  load 3 0.0 -2000.0 -168074.0
  load 4 0.0 -2000.0 168074.0
}
```

source example.tcl



Gravity Load Analysis

```
# first source in the model
source example.tcl

# create analysis & perform analysis
constraints transformation
numberer RCM
system BandGeneral
test NormDispIncr 1.0e-6 6 2
algorithm Newton
integrator LoadControl 0.1
analysis Static
analyze 10

# look at what happened to node 3
print node 3
```

* We will call this file example1.tcl for future examples

QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.

Cyclic Lateral Load Analysis

```
# first source in the model and do gravity load analysis  
source example1.tcl
```

```
# set gravity loads constant & reset time in domain  
loadConst -time 0.0
```

```
# create load pattern for lateral loads  
pattern Plain 2 Linear {  
    load 3 200.0 0.0 0.0  
    load 4 200.0 0.0 0.0  
}
```

```
# do some cyclic analysis  
foreach {numSteps stepSize} {10 0.1 10 -0.1 10 -0.1 10 0.1 10 0.1} {  
    integrator LoadControl $stepSize  
    analyze $numSteps  
  
    set time [getTime]  
    set disp [nodeDisp 3 1]  
    puts "Time: $time Displacement $disp"  
}
```

QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.

Transient Analysis - Uniform Excitation

```
# first source in the model and do gravity load analysis
source example1.tcl

# set gravity loads constant & reset time in domain
loadConst -time -0.01

# create load pattern
set motion IELC180
source READSMDFile.tcl
ReadSMDFile $motion.AT2 $motion.acc dt
set accelSeries "Path -filePath $motion.acc -dt $dt -factor 386.4"
pattern UniformExcitation 2 1 -accel $accelSeries

# set damping factors
rayleigh 0.0 0.0 0.0 0.0 0.0

# create the analysis
wipeAnalysis
System ProfileSPD
Test NormDispIncr 1.0e-16 10
Algorithm Newton
Integrator Newmark 0.5 0.25
Analysis Transient

# create a recorder
recorder Node -time -file example3.out -node 3 -dof 1 disp

analyze 4000 $dt
```

Transient Analysis - MultiSupport Excitation

```
# first source in the model and do gravity load analysis
source example1.tcl

# set gravity loads constant & reset time in domain
loadConst -time -0.01

# remove some sp constraints
remove sp 0
remove sp 3

# create load pattern
set motion IELC180
source READSMDFile.tcl
ReadSMDFile $motion.DT2 $motion.disp dt
set dispSeries "Path -filePath $motion.disp -dt $dt -factor 0.3937"
pattern MultiSupport 2 {
    groundMotion 5 Plain -disp $dispSeries
    imposedMotion 1 1 5
    imposedMotion 2 1 5
}


# set damping factors
rayleigh 0.0 0.0 0.0 0.0 0.0


# create the analysis
wipeAnalysis
system ProfileSPD
test NormDispIncr 1.0e-16 10
algorithm Newton
integrator Newmark 0.5 0.25
analysis Transient

analyze 4000 $dt

# set damping factors
rayleigh 0.0 0.0 0.0 0.0 0.0

# create a recorder
recorder Node -time -file example4.out -node 3 -dof 1 disp
```





Parameter Study - Response Spectra

```
source READSMDFile.tcl
modelBuilder BasicBuilder -ndm 1 -ndf 1
```

```
# set a bunch of parameters
set PI 3.14159265
set g 386.4
set TnMin 0.1; #min period
set TnMax 2.0; #max period
set TnIncr 0.1; #period incr
set M 1.0;      #mass
set A 1.0;      #area
set L 1.0;      #length
set motion ELCENTRO
set outFileID spectrum.dat
```

```
# open output file
Set outFileID [open $outFileID w]
```

```
#create accel series
ReadSMDFile $motion.AT2 $motion.acc dt
Set accelSeries "Path -filePath $motion.acc \
-dt $dt -factor $g"
```

```
# loop over period range
Set Tn $TnMin
while {$Tn <= $TnMax} {
  wipe
  set w [expr 2.0 * $PI / $Tn]
  set K [expr $w * $w * $M]
  set E [expr $k * $I / $A
```

```
node 1 0.0
node 2 $I -mass $M
fix 1 1
uniaxialMaterial Elastic 1 $E
element truss 1 1 2 $A 1
pattern UniformExcitation 2 1 -accel $accelSeries
rayleigh 0.0 0.0 0.0 0.0 0.0
```

```
recorder EnvelopeNode -file envelope.out -node 2 -dof 1 disp
system ProfileSPD
test NormDispIncr 1.0e-16 10
algorithm Newton
integrator Newmark 0.5 0.25
analysis Transient
analyze 2000 $dt
```

```
if [catch {open envelope.out r} inFileID]
  puts puts "ERROR - could not open file"
```

```
set min [gets $inFileID]
set max [gets $inFileID]
set absMax [gets $inFileID]
close $inFileID
puts outFileID "$Tn $absmax"
set Tn [expr $Tn + $TnIncr]
}
close $outFileID
```

QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.